

小鼠抗 TUBB3 单克隆抗体

中文名称: 小鼠抗 TUBB3 单克隆抗体

英文名称: Anti-TUBB3 mouse monoclonal antibody

别 名: tubulin beta 3 class III; CDCBM; FEOM3; TUBB4; CDCBM1; CFEOM3; beta-4;

CFEOM3A

相关类别: 内参抗体

储 存: 冷冻(-20℃)

抗 原: TUBB3

宿 主: Mouse

反应种属: Human, Mouse, Rat

标 记 物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:

This gene encodes a class III member of the beta tu bulin protein family. Beta tubulins are one of two co re protein families (alpha and beta tubulins) that het erodimerize and assemble to form microtubules. This protein is primarily expressed in neurons and may be involved in neurogenesis and axon guidance and maintenance. Mutations in this gene are the cause of congenital fibrosis of the extraocular muscles type 3. Alternate splicing results in multiple transcript variant s. A pseudogene of this gene is found on chromoso



	me 6.
Applications:	WB, IHC
Name of antibody:	TUBB3
Immunogen:	Fusion protein of human TUBB3
Full name:	tubulin beta 3 class III
Synonyms:	CDCBM; FEOM3; TUBB4; CDCBM1; CFEOM3; beta-4; C FEOM3A
SwissProt:	Q13509
IHC positive control:	Human prostate cancer and Human liver cancer; Hu man colorectal cancer
IHC Recommend dilution:	250-1000
WB Predicted band size:	50 kDa
WB Positive control:	Mouse brain tissue, Rat brain tissue and Human cere bella tissue lysates
WB Recommended dilution:	5000-20000







