

兔抗 RBFOX1 多克隆抗体

中文名称：兔抗 RBFOX1 多克隆抗体

英文名称： Anti-RBFOX1 rabbit polyclonal antibody

别名： 2BP1; FOX1; A2BP1; FOX-1; HRNBP1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RBFOX1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The Fox-1 family of RNA-binding proteins is evolutionarily conserved, and regulates tissue-specific alternative splicing in metazoa. Fox-1 recognizes a (U)GCAUG stretch in regulated exons or in flanking introns. The protein binds to the C-terminus of ataxin-2 and may contribute to the restricted pathology of spinocerebellar ataxia type 2 (SCA2). Ataxin-2 is the product of the SCA2 gene which causes familial neurodegenerative diseases. Fox-1 and ataxin-2 are both localized in the trans-Golgi network. Several alternatively spliced transcript variants encoding different isoforms have been found.

	nd for this gene.
Applications:	ELISA, IHC
Name of antibody:	RBFOX1
Immunogen:	Synthetic peptide of human RBFOX1
Full name:	RNA binding fox-1 homolog 1
Synonyms:	2BP1; FOX1; A2BP1; FOX-1; HRNBP1
SwissProt:	Q9NWB1
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	10-50

