

兔抗 NDUFS5 多克隆抗体

中文名称：兔抗 NDUFS5 多克隆抗体

英文名称： Anti-NDUFS5 rabbit polyclonal antibody

别名： CI15K; CI-15k

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NDUFS5

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of the NADH dehydrogenase (ubiquinone) iron-sulfur protein family. The encoded protein is a subunit of the NADH:ubiquinone oxidoreductase (complex I), the first enzyme complex in the electron transport chain located in the inner mitochondrial membrane. Alternative splicing results in multiple transcript variants and pseudogenes have been identified on chromosomes 1, 4 and 17.
Applications:	ELISA, WB, IHC
Name of antibody:	NDUFS5
Immunogen:	Fusion protein of human NDUFS5

Full name:	NADH dehydrogenase (ubiquinone) Fe-S protein 5, 15kDa (NADH-coenzyme Q reductase)
Synonyms :	CI15K; CI-15k
SwissProt:	O43920
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human brain
IHC Recommend dilution:	50-200
WB Predicted band size:	13 kDa
WB Positive control:	Human liver tissue and K562 cells, human fetal brain tissue and Raji cells
WB Recommended dilution:	500-2000



