

AP 标记的驴抗人 IgG

- 中文名称: AP标记的驴抗人 IgG
- 英文名称: AP-conjugated Donkey Anti-Human IgG
- 相关类别: AP标记二抗
- 储 存: 冷冻 (-20℃)
- 标记物: AP

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts wit h whole molecule human IgG. It also reacts with the light chains of o ther human immunoglobulins. No antibody was detected against non- immunoglobulin serum proteins. The antibody has been tested by ELI SA and/or solid-phase adsorbed to ensure minimal cross-reaction with bovine, chicken, goat, guinea pig, Syrian hamster, horse, mouse, rabbit , rat, and sheep serum proteins, but it may cross-react with immunogl obulins from other species.
Applications	1/2000-20000 for ELISA 1/2000-20000 for Western blotting