

兔抗 TUBGCP6 多克隆抗体

中文名称：兔抗 TUBGCP6 多克隆抗体

英文名称： Anti-TUBGCP6 rabbit polyclonal antibody

别名： GCP6; GCP-6; MCCR; MCPHCR

抗原： TUBGCP6

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene is part of a large multimeric complex required for microtubule nucleation at the centrosome.
Applications:	WB, IHC
Name of antibody:	TUBGCP6
Immunogen:	Synthesized peptide derived from internal of human TUBGCP6.
Full name:	tubulin, gamma complex associated protein 6
Synonyms :	GCP6; GCP-6; MCCR; MCPHCR
SwissProt:	Q96RT7
IHC positive control:	Human brain tissue

IHC Recommend dilution:	50-100
WB Predicted band size:	200 kDa
WB Positive control:	COLO cells lysate
WB Recommended dilution:	500-3000



