

AP 标记的驴抗山羊 IgG

中文名称：AP 标记的驴抗山羊 IgG

英文名称：AP-conjugated Donkey Anti-Goat IgG

相关类别：AP 标记二抗

储 存： 冷冻（-20℃）

标 记 物： AP

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule goat IgG. It also reacts with the light chains of other goat immunoglobulins. No antibody was detected against non - immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with chicken, guinea pig, Syrian hamster, horse, human, mouse, rabbit and rat serum proteins, but it may cross-react with immunoglobulins from other species.
Applications	1/2000-20000 for ELISA 1/2000-20000 for Western blotting