

AP 标记的驴抗兔 IgG

中文名称: AP 标记的驴抗兔 IgG

英文名称: AP-conjugated Donkey Anti-Rabbit IgG

相关类别: AP 标记二抗

储 存:冷冻(-20℃)

标 记 物: AP

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts w ith whole molecule rabbit IgG. It also reacts with the light chains of other rabbit immunoglobulins. No antibody was detected against no n- immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with bovine, chicken, goat, guinea pig, Syrian hamster, horse, human, mouse, rat, and sheep serum proteins, but it may cross-react with im munoglobulins from other species.
Applications	1/2000-20000 for ELISA 1/2000-20000 for Western blotting