

AP 标记的驴抗小鼠 IgG

- 中文名称: AP 标记的驴抗小鼠 IgG
- 英文名称: AP-conjugated Donkey Anti-Mouse IgG
- 相关类别: AP 标记二抗
- 储 存: 冷冻 (-20℃)
- 标记物: AP

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody react s with whole molecule mouse IgG. It also reacts with the light ch ains of other mouse immunoglobulins. No antibody was detected against non- immunoglobulin serum proteins. The antibody has b een tested by ELISA and/or solid-phase adsorbed to ensure mini mal cross-reaction with bovine, chicken, goat, guinea pig, Syrian h amster, horse, human, rabbit, rat, and sheep serum proteins, but i t may cross-react with immunoglobulins from other species.
Applications	1/2000-20000 for ELISA 1/2000-20000 for Western blotting