

兔抗 TUFM 多克隆抗体

中文名称：兔抗 TUFM 多克隆抗体

英文名称： Anti-TUFM rabbit polyclonal antibody

别名： P43; EFTU; COXPD4; EF-TuMT

抗原： TUFM

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a protein which participates in protein translation in mitochondria. Mutations in this gene have been associated with combined oxidative phosphorylation deficiency resulting in lactic acidosis and fatal encephalopathy. A pseudogene has been identified on chromosome 17.
Applications:	WB, IHC
Name of antibody:	TUFM
Immunogen:	Synthesized peptide derived from internal of human TUFM.
Full name:	Tu translation elongation factor, mitochondrial

Synonyms :	P43; EFTU; COXPD4; EF-TuMT
SwissProt:	P49411
IHC positive control:	Human lung carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	50 kDa
WB Positive control:	LOVO cells and NIH/3T3 cells lysates
WB Recommended dilution:	500-3000



