

## 兔抗 TULP2 多克隆抗体

中文名称：兔抗 TULP2 多克隆抗体

英文名称： Anti-TULP2 rabbit polyclonal antibody

别名： CT65; TUBL2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TULP2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	TULP2 is a member of a family of tubby-like genes (TULPs) that encode proteins of unknown function. Members of this family have been identified in plants, vertebrates, and invertebrates. The TULP proteins share a conserved C-terminal region of approximately 200 amino acid residues.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TULP2
<b>Immunogen:</b>	Fusion protein of human TULP2
<b>Full name:</b>	tubby like protein 2
<b>Synonyms :</b>	CT65; TUBL2

<b>SwissProt:</b>	O00295
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer and human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100

