

AP 标记山羊抗小鼠 IgG

中文名称: AP 标记山羊抗小鼠 IgG

英文名称: AP-conjugated Goat anti-mouse IgG

相关类别: AP 标记二抗

储 存: 冷冻(-20℃)

标 记 物: AP

技术规格

Background	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule mouse IgG. It also reacts with the light chains of other mouse immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.
Applications	1/2000-20000 for ELISA 1/2000-20000 for Western blotting