

兔抗 TUSC2 多克隆抗体

中文名称：兔抗 TUSC2 多克隆抗体

英文名称： Anti-TUSC2 rabbit polyclonal antibody

别名： tumor suppressor 2, mitochondrial calcium regulator; PAP; FUS1; PDAP2; C3orf11

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TUSC2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available.
Applications:	ELISA, IHC
Name of antibody:	TUSC2
Immunogen:	Fusion protein of human TUSC2
Full name:	tumor suppressor 2, mitochondrial calcium regulator
Synonyms:	PAP; FUS1; PDAP2; C3orf11
SwissProt:	O75896
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human lung cancer
IHC Recommend dilution:	200-400

