

## AP 标记兔抗人 IgG

中文名称: AP 标记兔抗人 IgG

英文名称: AP-conjugated Rabbit anti-human IgG

相关类别: AP 标记二抗

储 存: 冷冻 (-20℃)

标 记 物: AP

### 技术规格

<b>Background</b>	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule human IgG. It also reacts with the light chains of other human immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.
<b>Applications</b>	1/2000-20000 for ELISA 1/2000-20000 for Western blotting