

## GPBP1L1 抗原（重组蛋白）

中文名称：GPBP1L1 抗原（重组蛋白）

英文名称：GPBP1L1 Antigen (Recombinant Protein)

别名：GC-rich promoter binding protein 1 like 1; SP192

储存：冷冻（-20℃）

相关类别：抗原

### 概述

Fusion protein corresponding to a region derived from 275-474 amino acids of human GPBP1L1

### 技术规格

|                           |  |
|---------------------------|--|
| <b>Full name:</b>         | GC-rich promoter binding protein 1 like 1              |
| <b>Synonyms:</b>          | SP192  |
| <b>Swissprot:</b>         | Q9HC44   |
| <b>Gene Accession:</b>    | BC015713   |
| <b>Purity:</b>            | >85%, as determined by Coomassie blue stained SDS-PAGE |
| <b>Expression system:</b> | Escherichia coli                                       |
| <b>Tags:</b>              | His tag C-Terminus, GST tag N-Terminus                 |
| <b>Background:</b>        | Possible transcription factor.                         |