

## 兔抗 TWF1 多克隆抗体

中文名称：兔抗 TWF1 多克隆抗体

英文名称： Anti-TWF1 rabbit polyclonal antibody

别名： A6; PTK9

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TWF1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes twinfilin, an actin monomer-binding protein conserved from yeast to mammals. Studies of the mouse counterpart suggest that this protein may be an actin monomer-binding protein, and its localization to cortical G-actin-rich structures may be regulated by the small GTPase RAC1.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	TWF1
<b>Immunogen:</b>	Fusion protein of human TWF1
<b>Full name:</b>	twinfilin actin-binding protein 1
<b>Synonyms :</b>	A6; PTK9

<b>SwissProt:</b>	Q12792
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human esophagus cancer and human brain
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	40 kDa
<b>WB Positive control:</b>	Lovo cells
<b>WB Recommended dilution:</b>	200-1000



