

兔抗 NECTIN2 多克隆抗体

中文名称：兔抗 NECTIN2 多克隆抗体

英文名称：Anti-NECTIN2 rabbit polyclonal antibody

别名：nectin cell adhesion molecule 2; HVEB; PRR2; CD112

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NECTIN2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a single-pass type I membrane glycoprotein with two Ig-like C2-type domains and an Ig-like V-type domain. This protein is one of the plasma membrane components of adherens junctions. It also serves as an entry for certain mutant strains of herpes simplex virus and pseudorabies virus, and it is involved in cell to cell spreading of these viruses. Variations in this gene have been associated with differences in the severity of multiple sclerosis. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]
Applications:	ELISA, IHC

Name of antibody:	NECTIN2
Immunogen:	Fusion protein of human NECTIN2
Full name:	nectin cell adhesion molecule 2
Synonyms:	HVEB; PRR2; CD112; PVRL2; PVRR2
SwissProt:	Q92692
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200

