

兔抗 RBM17 多克隆抗体

中文名称：兔抗 RBM17 多克隆抗体

英文名称：Anti-RBM17 rabbit polyclonal antibody

别名：RNA binding motif protein 17; SPF45

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RBM17

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an RNA binding protein. The encoded protein is part of the spliceosome complex and functions in the second catalytic step of mRNA splicing. Alternatively spliced transcript variants have been described. Related pseudogenes exist on chromosomes 9 and 15.
Applications:	ELISA, WB, IHC
Name of antibody:	RBM17
Immunogen:	Fusion protein of human RBM17
Full name:	RNA binding motif protein 17
Synonyms:	SPF45

SwissProt:	Q96I25
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	45 kDa
WB Positive control:	Mouse liver tissue lysate
WB Recommended dilution:	200-1000



