

兔抗 IFI27 多克隆抗体

中文名称：兔抗 IFI27 多克隆抗体

英文名称：Anti-IFI27 rabbit polyclonal antibody

别名：interferon alpha inducible protein 27; P27; ISG12; FAM14D; ISG12A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IFI27

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Promotes cell death. Mediates IFN-induced apoptosis characterized by a rapid and robust release of cytochrome C from the mitochondria and activation of BAX and caspases 2, 3, 6, 8 and 9.
Applications:	ELISA, IHC
Name of antibody:	IFI27
Immunogen:	Fusion protein of human IFI27
Full name:	interferon alpha inducible protein 27
Synonyms:	P27; ISG12; FAM14D; ISG12A
SwissProt:	P40305
ELISA Recommended dilu	5000-10000

tion:	
IHC positive control:	Human tonsil
IHC Recommend dilution:	100-200

