

兔抗 TWF2 多克隆抗体

中文名称：兔抗 TWF2 多克隆抗体

英文名称： Anti-TWF2 rabbit polyclonal antibody

别名： A6r; A6RP; PTK9L; MSTP011

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： TWF2

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene was identified by its interaction with the catalytic domain of protein kinase C-zeta. The encoded protein contains an actin-binding site and an ATP-binding site. It is most closely related to twinfilin (PTK9), a conserved actin monomer-binding protein.
Applications:	ELISA, IHC
Name of antibody:	TWF2
Immunogen:	Fusion protein of human TWF2
Full name:	twinfilin actin-binding protein 2
Synonyms :	A6r; A6RP; PTK9L; MSTP011

SwissProt:	Q6IBS0
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

