

GST 抗原（重组蛋白）

中文名称： GST 抗原（重组蛋白）

英文名称： GST Antigen (Recombinant Protein)

相关类别： 抗原

储 存： 超低温（-80℃）

概述

Full length protein of GST

技术规格

Full name:	Glutathione S-Transferase
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Background:	GST (Glutathione S-Transferase) is a 26kDa protein which was widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins. GST isoenzymes appear to play a central role in the parasite detoxification system. Other functions are also suspected including a role in increasing the solubility of haematin in the parasite gut.