

兔抗 TWIST1 多克隆抗体

中文名称：兔抗 TWIST1 多克隆抗体

英文名称： Anti-TWIST1 rabbit polyclonal antibody

别名： SCS; ACS3; CRS1; BPES2; BPES3; TWIST; bHLHa38

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TWIST1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Basic helix-loop-helix (bHLH) transcription factors have been implicated in cell lineage determination and differentiation. The protein encoded by this gene is a bHLH transcription factor and shares similarity with another bHLH transcription factor, Dermo1. The strongest expression of this mRNA is in placental tissue; in adults, mesodermally derived tissues express this mRNA preferentially. Mutations in this gene have been found in patients with Saethre-Chotzen syndrome.
Applications:	ELISA, IHC

Name of antibody:	TWIST1
Immunogen:	Synthetic peptide of human TWIST1
Full name:	twist basic helix-loop-helix transcription factor 1
Synonyms :	SCS; ACS3; CRS1; BPES2; BPES3; TWIST; bHLHa38
SwissProt:	Q15672
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and Human liver cancer
IHC Recommend dilution:	25-100

