

兔抗 RBM26 多克隆抗体

中文名称：兔抗 RBM26 多克隆抗体

英文名称：Anti-RBM26 rabbit polyclonal antibody

别名：RNA binding motif protein 26; ARRS2; SE70-2; ZC3H17; PRO1777; C13orf10; PPP1R132

相关类别：一抗

抗原：RBM26

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	RBM26 (RNA binding motif protein 26) bears two RNA recognition motif (RRM) domains and a C3H1-type zinc finger. The function of RBM26 has not been characterized.
Applications:	ELISA, IHC
Name of antibody:	RBM26
Immunogen:	Synthetic peptide of human RBM26
Full name:	RNA binding motif protein 26
Synonyms:	ARRS2; SE70-2; ZC3H17; PRO1777; C13orf10; PPP1R132

SwissProt:	Q5T8P6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	20-100

