

GTF2F1 抗原（重组蛋白）

中文名称：GTF2F1 抗原（重组蛋白）

英文名称：GTF2F1 Antigen (Recombinant Protein)

别名：general transcription factor IIF subunit 1; BTF4; RAP74; TF2F1; TFIIF

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 318-517 amino acids of human GTF2F1

技术规格

Full name:	general transcription factor IIF subunit 1
Synonyms:	BTF4; RAP74; TF2F1; TFIIF
Swissprot:	P35269
Gene Accession:	BC013007
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	TFIIF is a general transcription initiation factor that binds to RNA polymerase II and helps to recruit it to the initiation complex in collaboration with TFIIB. It promotes transcription elongation.