

## 兔抗 RBMX 多克隆抗体

中文名称：兔抗 RBMX 多克隆抗体

英文名称： Anti-RBMX rabbit polyclonal antibody

别名： RNMX; HNRPG; HNRNPG; RBMXP1; RBMXRT; hnRNP-G

相关类别： 一抗

抗原： RBMX

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

This gene belongs to the RBMY gene family which includes candidate Y chromosome spermatogenesis genes. This gene, an active X chromosome homolog of the Y chromosome RBMY gene, is widely expressed whereas the RBMY gene evolved a male-specific function in spermatogenesis. Pseudogenes of this gene, found on chromosomes 1, 4, 9, 11, and 6, were likely derived by retrotransposition from the original gene. Alternatively spliced transcript variants encoding different isoforms have been identified. A snoRNA gene (SNORD61) is found in one of its introns.

<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	RBMX
<b>Immunogen:</b>	Fusion protein of human RBMX
<b>Full name:</b>	RNA binding motif protein, X-linked
<b>Synonyms :</b>	RNMX; HNRPG; HNRNPG; RBMXP1; RBMXRT; hnRNP-G
<b>SwissProt:</b>	P38159
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	42 kDa
<b>WB Positive control:</b>	K562, hela and Jurkat cells
<b>WB Recommended dilution:</b>	200-1000

