

兔抗 NEGR1 多克隆抗体

中文名称：兔抗 NEGR1 多克隆抗体

英文名称：Anti-NEGR1 rabbit polyclonal antibody

别名：neuronal growth regulator 1; Ntra; KILON; IGLON4; DMML2433

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NEGR1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in cell-adhesion. May function as a trans-neural growth-promoting factor in regenerative axon sprouting in the mammalian brain (By similarity).
Applications:	ELISA, IHC
Name of antibody:	NEGR1
Immunogen:	Fusion protein of human NEGR1
Full name:	neuronal growth regulator 1
Synonyms:	Ntra; KILON; IGLON4; DMML2433
SwissProt:	Q7Z3B1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil

IHC Recommend dilution:	50-300
-------------------------	--------

