

## 兔抗 NEIL1 多克隆抗体

中文名称：兔抗 NEIL1 多克隆抗体

英文名称： Anti-NEIL1 rabbit polyclonal antibody

别名： FPG1; NEI1; hFPG1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NEIL1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the Nei endonuclease VIII-like gene family which encodes DNA glycosylases. The encoded enzyme participates in the DNA repair pathway by initiating base excision repair by removing damaged bases, primarily oxidized pyrimidines. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NEIL1
<b>Immunogen:</b>	Fusion protein of human NEIL1
<b>Full name:</b>	nei like DNA glycosylase 1

<b>Synonyms:</b>	FPG1; NEI1; hFPG1
<b>SwissProt:</b>	Q96FI4
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer;Human brain
<b>IHC Recommend dilution:</b>	25-100



