

兔抗 RBP4 多克隆抗体

中文名称：兔抗 RBP4 多克隆抗体

英文名称： Anti-RBP4 rabbit polyclonal antibody

别名： RDCCAS

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RBP4

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This protein belongs to the lipocalin family and is the specific carrier for retinol (vitamin A alcohol) in the blood . It delivers retinol from the liver stores to the peripheral tissues. In plasma, the RBP-retinol complex interacts with transthyretin which prevents its loss by filtration through the kidney glomeruli. A deficiency of vitamin A blocks secretion of the binding protein posttranslationally and results in defective delivery and supply to the epidermal cells.
Applications:	ELISA, WB
Name of antibody:	RBP4

Immunogen:	Fusion protein of human RBP4
Full name:	retinol binding protein 4, plasma
Synonyms :	RDCCAS
SwissProt:	P02753
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	23 kDa
WB Positive control:	Human fetal liver and liver cancer tissue
WB Recommended dilution:	200-1000

