

## 生物素标记驴抗人 IgG

中文名称: 生物素标记驴抗人 IgG

英文名称: Biotin-conjugated Donkey Anti-Human IgG

相关类别: Biotin 标记二抗

储 存: 冷冻(-20℃)

标 记 物: Biotin

## 技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with h whole molecule human IgG. It also reacts with the light chains of o ther human immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody has been tested by ELI SA and/or solid-phase adsorbed to ensure minimal cross-reaction with bovine, chicken, goat, guinea pig, Syrian hamster, horse, mouse, rabb it, rat, and sheep serum proteins, but it may cross-react with immuno globulins from other species.
Applications	1/10000-200000 for ELISA 1/200-2000 for enzyme immunohisto/cytoch emistry 1/100-500 for Flow Cytometer and fluorescence immunohisto/cytochemistry