

生物素标记的驴抗小鼠 IgG

中文名称： 生物素标记的驴抗小鼠 IgG

英文名称： Biotin-conjugated Donkey Anti-Mouse IgG

相关类别： Biotin 标记二抗

储 存： 冷冻（-20℃）

标 记 物： Biotin

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule mouse IgG. It also reacts with the light chains of other mouse immunoglobulins. No antibody was detected against non- immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with bovine, chicken, goat, guinea pig, Syrian hamster, horse, human, rabbit, rat, and sheep serum proteins, but it may cross-react with immunoglobulins from other species.
Applications	1/10000-200000 for ELISA 1/200-2000 for enzyme immunohisto/chemistry 1/100-500 for Flow Cytometer and fluorescence immunohisto/chemistry