

## 兔抗 TXN2 多克隆抗体

中文名称：兔抗 TXN2 多克隆抗体

英文名称： Anti-TXN2 rabbit polyclonal antibody

别名： MTRX; TRX2; MT-TRX

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TXN2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This nuclear gene encodes a mitochondrial member of the thioredoxin family, a group of small multifunctional redox-active proteins. The encoded protein may play important roles in the regulation of the mitochondrial membrane potential and in protection against oxidant-induced apoptosis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TXN2
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	thioredoxin 2
<b>Synonyms :</b>	MTRX; TRX2; MT-TRX
<b>SwissProt:</b>	Q99757

<b>ELISA Recommended dilution</b> :	2000-5000
<b>IHC positive control:</b>	Human liver cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100

