

HCRTR2 抗原（多肽片段）

中文名称： HCRTR2 抗原（多肽片段）

英文名称： HCRTR2 Antigen

相关类别： 抗原

储 存： 冷冻（-20℃）

概述

Synthetic peptide corresponding to a region derived from 405-419 amino acids (CKSLTTQISNFDNISK) of Human HCRTR2

技术规格

Full name:	hypocretin (orexin) receptor 2
Synonyms:	OX2R
Gene Accession:	NP_001517
Swissprot:	O43614
Background:	The protein encoded by this gene is a G-protein coupled receptor involved in the regulation of feeding behavior. The encoded protein binds the hypothalamic neuropeptides orexin A and orexin B. A related gene (HCRTR1) encodes a G-protein coupled receptor that selectively binds orexin A.