

生物素标记山羊抗人 IgM

中文名称: 生物素标记山羊抗人 IgM

英文名称: Biotin-conjugated Goat anti-human IgM

相关类别: Biotin 标记二抗

储 存: 冷冻(-20℃)

技术规格

Background	IgM normally constitutes about 10% of serum immunoglobulins. I gM antibody is prominent in early immune responses to most an tigens and is largely confined to plasma due to its large size. Mo nomeric IgM is expressed as a membrane bound antibody on the surface of B cells and as a pentamer when secreted by plasma ce lls. Due to its high valency IgM is more efficient than other isoty pes is binding antigens with repeating epitopes (virus particles an d red blood cells) and is more efficient than IgG in activiating the complement pathway. The gene for the mu constant region contains four domains separated by short intervening sequences.
Applications	Optimal tested applications and working dilutions must be deter mined by the end user.