

Hela 细胞裂解液

中文名称： Hela 细胞裂解液

英文名称： Hela Cell Lysate

别名： Hela 细胞裂解液

相关类别： 抗原

储存： 冷冻（-20℃） 避光

概述

Hela cell lysates from human cervical cancer cells were cultured in the media containing 10% FBS.

技术规格

Backg round	Cells were lysed in lysis buffer (50mM Tris(PH7.4), 150mM NaCl, 1% NP-40, 1% SDS). 100 µg of cellular protein was mixed with 4 times volume of 5X SDS-PAGE loading buffer.
------------------------	--