

兔抗 NEK7 多克隆抗体

中文名称: 兔抗 NEK7 多克隆抗体

英文名称: Anti-NEK7 rabbit polyclonal antibody

抗 原: NEK7

储 存: 冷冻 (-20℃)

宿 主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	NIMA-related kinases share high amino acid sequence identity with the gene product of the <i>Aspergillus nidulans</i> 'n ever in mitosis A' gene, which controls initiation of mitosis. Protein kinase which plays an important role in mitotic cell cycle progression. Required for microtubule nucleation activity of the centrosome, robust mitotic spindle formation and cytokinesis. Phosphorylates RPS6KB1.
Applications:	ELISA, IHC
Name of antibody:	NEK7
Immunogen:	Synthetic peptide of human NEK7
Full name:	NIMA related kinase 7
SwissProt:	Q8TDX7
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human esophagus cancer and human breast cancer
IHC Recommend dilution:	50-300



