

## 兔抗 IFNA2 多克隆抗体

中文名称: 兔抗 IFNA2 多克隆抗体

英文名称: Anti-IFNA2 rabbit polyclonal antibody

别名: interferon alpha 2; IFNA; INFA2; IFNA2B; IFN-alphaA

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: IFNA2

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the alpha interferon gene cluster on chromosome 9. The encoded protein is a cytokine produced in response to viral infection. Use of the recombinant form of this protein has been shown to be effective in reducing the symptoms and duration of the common cold.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	IFNA2
<b>Immunogen:</b>	Synthetic peptide of human IFNA2
<b>Full name:</b>	interferon alpha 2
<b>Synonyms:</b>	IFNA; INFA2; IFNA2B; IFN-alphaA
<b>SwissProt:</b>	P01563

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	30-150



