

HEPG2.2.15 细胞裂解液

中文名称： HEPG2.2.15 细胞裂解液

英文名称： HEPG2.2.15 Cell Lysate

别名： HEPG2.2.15 细胞裂解液

相关类别： 抗原

储存： 冷冻（-20℃） 避光

概述

HEPG2.2.15 cell lysates from human hepatoma cells were cultured in the media containing 20% FBS.

技术规格

Background	Cells were lysed in lysis buffer (50mM Tris(PH7.4), 150mM NaCl, 1% NP-40, 1% SDS). 100 µg of cellular protein was mixed with 4 times volume of 5X SDS-PAGE loading buffer.
-------------------	--