

生物素标记小鼠抗兔 IgG

中文名称： 生物素标记小鼠抗兔 IgG

英文名称： Biotin-conjugated Mouse anti-rabbit IgG

相关类别： Biotin 标记二抗

储 存： 冷冻（-20℃）

标记物： Biotin

技术规格

Background	<p>Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule rabbit IgG. It also reacts with the light chains of other rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins.</p> <p>The antibody exhibits inherent minimal cross-reaction to mouse serum proteins and has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with human, goat, mouse, and sheep serum proteins. The antibody may cross-react with immunoglobulins from other species.</p>
Applications	1/5000-100000 for ELISA and Western blotting using enzyme-conj

ugated streptavidin 1/100-1000 for enzyme immunohisto/cytoche
mistry 1/50-200 for Flow Cytometer and fluorescence immunohist
o/cytochemistry