

兔抗 NEK8 多克隆抗体

中文名称：兔抗 NEK8 多克隆抗体

英文名称： Anti-NEK8 rabbit polyclonal antibody

别名： JCK; NPHP9; RHPD2; NEK12A

相关类别： 一抗

抗原： NEK8

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the serine/threonine protein kinase family related to NIMA (never in mitosis, gene A) of <i>Aspergillus nidulans</i> . The encoded protein may play a role in cell cycle progression from G2 to M phase. Mutations in the related mouse gene are associated with a disease phenotype that closely parallels the juvenile autosomal recessive form of polycystic kidney disease in humans.
Applications:	ELISA, IHC
Name of antibody:	NEK8
Immunogen:	Fusion protein of human NEK8

Full name:	NIMA-related kinase 8
Synonyms :	JCK; NPHP9; RHPD2; NEK12A
SwissProt:	Q86SG6
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	25-100

