

HRP 标记的驴抗鸡 IgG

中文名称： HRP 标记的驴抗鸡 IgG

英文名称： HRP-conjugated Donkey Anti-Chicken IgG

相关类别： HRP 标记二抗

储 存： 冷冻 (-20℃)

标 记 物： HRP

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule Chicken IgY. It also reacts with the light chains of other Chicken immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with bovine, goat, guinea pig, Syrian hamster, horse, human, mouse, rabbit, rat, and sheep serum proteins, but it may cross-react with immunoglobulins from other species.
Applications	1/2000-50000 for ELISA 1/200-2000 for Immunohisto/cytochemistry 1/5000-100000 for Western blotting with ECL substrates