

## 兔抗 NEK9 多克隆抗体

中文名称：兔抗 NEK9 多克隆抗体

英文名称：Anti-NEK9 rabbit polyclonal antibody

别名：Nek8; NERCC; NERCC1

相关类别：一抗 储存 冷冻 (-20℃)

宿主：Rabbit

抗原：NEK9

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	NEK9, a NEK type protein kinase, regulates chromosome alignment and segregation in mitosis. The protein has a N-terminal NIMA-like catalytic domain, a central domain with homology to the guanine nucleotide exchange factor for the GTPase Ran (RCC1), and a C-terminal coiled-coil domain. It is phosphorylated by active p34(Cdc2) and is capable of autophosphorylation and oligomerization.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NEK9
<b>Immunogen:</b>	Synthetic peptide of human NEK9
<b>Full name:</b>	NIMA-related kinase 9
<b>Synonyms :</b>	Nek8; NERCC; NERCC1
<b>SwissProt:</b>	Q8TD19

<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human cervical cancer and Human breast cancer
<b>IHC Recommend dilution:</b>	25-100

