

兔抗 IFNA8 多克隆抗体

中文名称：兔抗 IFNA8 多克隆抗体

英文名称：Anti-IFNA8 rabbit polyclonal antibody

别名：interferon alpha 8; IFN-alphaB

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IFNA8

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.
Applications:	ELISA, WB
Name of antibody:	IFNA8
Immunogen:	Synthetic peptide of human IFNA8
Full name:	interferon alpha 8
Synonyms:	IFN-alphaB
SwissProt:	P32881
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	22 kDa

WB Positive control:	Human fetal liver tissue lysate
WB Recommended dilution:	500-2000

