

## 兔抗 IFNA16 多克隆抗体

中文名称：兔抗 IFNA16 多克隆抗体

英文名称：Anti-IFNA16 rabbit polyclonal antibody

别名：IFN-alphaO

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IFNA16

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Interferon, alpha 16 is produced by macrophages and has antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	IFNA16
<b>Immunogen:</b>	Synthetic peptide of human IFNA16
<b>Full name:</b>	interferon, alpha 16
<b>Synonyms :</b>	IFN-alphaO
<b>SwissProt:</b>	P05015
<b>ELISA Recommended dilution:</b>	2000-5000

<b>WB Predicted band size:</b>	22 kDa
<b>WB Positive control:</b>	Riji cells
<b>WB Recommended dilution:</b>	500-2000

