

HRP 标记的驴抗山羊 IgG

中文名称: HRP标记的驴抗山羊 IgG

英文名称: HRP-conjugated Donkey Anti-Goat IgG

相关类别: HRP 标记二抗

储存: 冷冻(-20℃) 避光

标记物: HRP

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody r eacts with whole molecule goat IgG. It also reacts with the lig ht chains of other goat immunoglobulins. No antibody was de tected against non- immunoglobulin serum proteins. The antib ody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with chicken, guinea pig, Syrian hamster, horse, human, mouse, rabbit and rat serum proteins, but it may cross-react with immunoglobulins from other species.
Applications	1/2000-50000 for ELISA 1/200-2000 for Immunohisto/cytochem istry 1/5000-100000 for Western blotting with ECL substrates