

HRP 标记的驴抗兔 IgG

中文名称: HRP 标记的驴抗兔 IgG

英文名称: HRP-conjugated Donkey Anti-Rabbit IgG

相关类别: HRP 标记二抗

储 存: 冷冻(-20℃) 避光

标 记 物: HRP

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule rabbit IgG. It also reacts with the light chains of other r abbit immunoglobulins. No antibody was detected against non- immuno globulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with bovine, chic ken, goat, guinea pig, Syrian hamster, horse, human, mouse, rat, and sh eep serum proteins, but it may cross-react with immunoglobulins from other species.
Applications	1/2000-50000 for ELISA 1/200-2000 for Immunohisto/cytochemistry 1/50 00-100000 for Western blotting with ECL substrates