

HRP 标记的驴抗豚鼠 lgG

- 中文名称: HRP标记的驴抗豚鼠 IgG
- 英文名称: HRP-conjugated Donkey Anti-Guinea pig IgG
- 相关类别: HRP标记二抗
- 储 存: 冷冻 (-20℃)
- 标 记 物: HRP

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule pig IgG. It also reacts with the light chains of other pig immunoglobulins. No antibody was detected against non- immunoglobuli n serum proteins. The antibody has been tested by ELISA and/or solid-p hase adsorbed to ensure minimal cross-reaction with bovine, chicken, go at, horse, human, mouse, rabbit, rat, and sheep serum proteins, but it m ay cross-react with immunoglobulins from other species.
Applications	1/2000-50000 for ELISA 1/200-2000 for Immunohisto/cytochemistry 1/500 0-100000 for Western blotting with ECL substrates