

## 兔抗 RDH11 多克隆抗体

中文名称：兔抗 RDH11 多克隆抗体

英文名称： Anti-RDH11 rabbit polyclonal antibody

别名： MDT1; CGI82; PSDR1; RALR1; SCALD; ARSDR1; HCBP12; SDR7C1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RDH11

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is an NADPH-dependent retinal reductase and a short-chain dehydrogenase/reductase. The encoded protein has no steroid dehydrogenase activity.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RDH11
<b>Immunogen:</b>	Fusion protein of human RDH11
<b>Full name:</b>	retinol dehydrogenase 11 (all-trans/9-cis/11-cis)
<b>Synonyms :</b>	MDT1; CGI82; PSDR1; RALR1; SCALD; ARSDR1; HCBP12; SDR7C1
<b>SwissProt:</b>	Q8TC12

<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human prostate cancer
<b>IHC Recommend dilution:</b>	10-50

