

兔抗 NET1 多克隆抗体

中文名称：兔抗 NET1 多克隆抗体

英文名称：Anti-NET1 rabbit polyclonal antibody

别名：NET1; ARHGEF8; NET1A

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：NET1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is part of the family of Rho guanine nucleotide exchange factors. Members of this family activate Rho proteins by catalyzing the exchange of GDP for GTP. The protein encoded by this gene interacts with RhoA within the cell nucleus and may play a role in repairing DNA damage after ionizing radiation. Pseudogenes of this gene are located on the long arms of chromosomes 1, 7 and 18. Alternative splicing results in multiple transcript variants that encode different protein isoforms.
Applications:	WB
Name of antibody:	NET1

Immunogen:	Fusion protein of human NET1
Full name:	neuroepithelial cell transforming 1
Synonyms :	NET1; ARHGEF8; NET1A
SwissProt:	Q7Z628
WB Predicted band size:	68 kDa
WB Positive control:	Jurkat and 293T cells
WB Recommended dilution:	500-2000

