

## 兔抗 TYRP1 多克隆抗体

- 中文名称：兔抗 TYRP1 多克隆抗体
- 英文名称：Anti-TYRP1 rabbit polyclonal antibody
- 别名：TRP; CAS2; CATB; GP75; OCA3; TRP1; TYRP; b-PROTEIN
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：TYRP1
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a melanosomal enzyme that belongs to the tyrosinase family and plays an important role in the melanin biosynthetic pathway. Defects in this gene are the cause of rufous oculocutaneous albinism and oculocutaneous albinism type III.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TYRP1
<b>Immunogen:</b>	Synthetic peptide of human TYRP1
<b>Full name:</b>	tyrosinase-related protein 1
<b>Synonyms:</b>	TRP; CAS2; CATB; GP75; OCA3; TRP1; TYRP; b-PROTEIN

<b>SwissProt:</b>	P17643
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer and human lung cancer
<b>IHC Recommend dilution:</b>	20-100

