

## 兔抗 IFNG 多克隆抗体

中文名称:	兔抗 IFNG 多克隆抗体
英文名称:	Anti-IFNG rabbit polyclonal antibody
别名:	IFG; IFI
抗原:	IFNG
储存:	冷冻 (-20℃)
宿主:	Rabbit
相关类别:	一抗
反应种属:	Human
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the type II interferon family. The protein encoded is a soluble cytokine with antiviral, immunoregulatory and anti-tumor properties and is a potent activator of macrophages. Mutations in this gene are associated with aplastic anemia.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	IFNG
<b>Immunogen:</b>	Fusion protein of human IFNG
<b>Full name:</b>	interferon, gamma
<b>Synonyms :</b>	IFG; IFI
<b>SwissProt:</b>	P01579

<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human cervical cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	19 kDa
<b>WB Positive control:</b>	Human fetal liver tissue
<b>WB Recommended dilution:</b>	200-1000





