

兔抗 REC8 多克隆抗体

中文名称：兔抗 REC8 多克隆抗体

英文名称： Anti-REC8 rabbit polyclonal antibody

别名： Rec8p; REC8L1; HR21spB

抗原： REC8

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the kleisin family of SMC (structural maintenance of chromosome) protein partners. The protein localizes to the axial elements of chromosomes during meiosis in both oocytes and spermatocytes. In the mouse, the homologous protein is a key component of the meiotic cohesion complex, which regulates sister chromatid cohesion and recombination between homologous chromosomes.
Applications:	ELISA, IHC
Name of antibody:	REC8

Immunogen:	Fusion protein of human REC8
Full name:	REC8 meiotic recombination protein
Synonyms :	Rec8p; REC8L1; HR21spB
SwissProt:	O95072
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human prostate cancer
IHC Recommend dilution:	20-100

